




Sealing Tape Pads for use with Empore™ 96 Well Plates

Suggested Instructions for Use

CDS Sealing Tape Pads can be used to cover unused wells during an extraction assay. Use a clean sheet each time for maximum sealing ability; these sheets are not intended for reuse.

NOTE: CDS Does not recommend the use of these sheets for sealing collection plates prior to sample injection because adhesive may interfere with analysis.



Customer Contact Information

Global Customer and Technical Service

Phone 1-800-541-6593

Fax 1-610-932-4158

Website www.CDSAnalytical.com/Empore



34-8711-5024-8_typed.indd 1

1/2/2013 8:50:33 AM



WARNING: To reduce the risks associated with improper disposal and/or handling of contaminants in used cards: Take appropriate steps to assess the disposal required for any altered product or materials added to the product. (Alteration of the product or the addition of other materials to the product may require different disposal methods.)

Warranty Information

All statements, technical information and recommendations herein are based on tests CDS believes to be reliable, but the accuracy or completeness thereof is not guaranteed. CDS WARRANTS ONLY THAT CDS PRODUCTS WILL MEET CDS'S SPECIFICATIONS AT THE TIME OF SHIPMENT. THE FOREGOING WARRANTY IS MADE IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND FREEDOM FROM NON-INFRINGEMENT. CDS's Empore™ Sample Preparation Products are intended for solid phase extraction during scientific research only. These products are not intended or warranted for use in medical devices or in assessment and treatment of clinical patients.

Limitations of Liabilities:

THE REMEDIES PROVIDED HEREIN ARE EXCLUSIVE REMEDIES AGAINST CDS FOR ANY ALLEGED OR ACTUAL NONCONFORMANCE TO SPECIFICATIONS OR DEFECT OR OTHER FAILURE IN PRODUCTS. UNDER NO CIRCUMSTANCES IS CDS LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES (INCLUDING LOST PROFITS) IN ANY WAY RELATED TO THE PRODUCT OR SUPPLY OF PRODUCT UNDER ANY THEORY OF LAW INCLUDING, BUT NOT LIMITED TO CONTRACT, NEGLIGENCE AND STRICT LIABILITY.



ABQ ?I_jwra_j* JJA
465 Limestone Road
Oxford, PA 19363-0277
Phone: 1-800-541-6593
Fax: 1-610-932-4158
Website: www.CDSAnalytical.com/Empore

CDS and Empore are trademarks of CDS.
34-8711-5024-8

34-8711-5024-8_typed.indd 2

1/2/2013 8:51:15 AM